



# TOUCH<sup>®</sup>

Dual mobility trapeziometacarpal prosthesis



# TOUCH®

## Dual mobility trapeziometacarpal prosthesis Class III



### CUPS


2 cup sizes (diameter 9 mm and 10 mm)

5 anti-rotation fins

Spiked crown for stability

Material: Stainless steel

Dual coating: Porous titanium + HA



### INSERTS

Only one pre-assembled insert size on each neck

Dual mobility retentive insert

Internal diameter 4 mm, external diameter 7 mm

Material: Polyethylene


### NECKS

3 heights: 6, 8 and 10 mm

2 types of neck: straight and 15° offset

Head diameter: 4 mm

Material: Stainless steel



### STEMS

5 stem models

4 diameters: 7.5, 8.5, 9.5 and 10.5 mm

Anatomical trilobe section

Cone for neck fixation

Anti-rotation neck system


Height: 25 mm for T0, 28 mm from T1 to T4

Material: Titanium

Dual coating: Porous titanium + HA



### INSTRUMENT SET



Supplied by:  
**Keri Medical Switzerland - HQ**  
Route des Acacias, 45a  
1227 Geneva - Switzerland  
Tel: +41 58 255 01 30  
Fax: +41 58 255 01 40  
customerservice@kerimedical.com

**Keri Medical France**  
ActiPro - 219 rue Laszlo Biro  
74160 Archamps - France  
Tel: +33(0)4 50 87 66 49  
Fax: +33 (0)4 80 16 05 10  
adv@kerimedical.com

[www.kerimedical.com](http://www.kerimedical.com)

