

IMPLATE[®]


wrist arthrodesis nail



Put optimal positioning in the patient's hands

Low-profile intramedullary construct tailored to the functional needs of the patient

Intramedullary placement minimizes soft tissue irritation




IMPLATE Wrist Fusion Nail



Alternate Wrist Fusion Plating Constructs



Modular design allows wrist positioning for power or digital dexterity



MicroSplines allow for 360° rotational adjustment



Threaded slots provide up to 8mm of translation for compression or distraction

3 METACARPAL NAILS
4.0mm Mini, 4.0mm, 4.6mm

2 RADIAL NAILS
Long, Short



PUBLISHED CLINICAL DATA:

Orbay JL, Feliciano E, Orbay C. (2012) Locked Intramedullary Total Wrist Arthrodesis. *J Wrist Surg*, 1(2): 179-184. DOI: 10.1055/s-0032-1329630