

— DISTAL RADIUS PROSTHESIS

COBRA™



DOCUMENTATION



groupe lépine

Depuis 1714

DISTAL RADIUS PROSTHESIS

- cobalt–chromium–molybdenum alloy.
- cementless fixation made possible by bilayer coating of vacuum-plasma sprayed pure porous titanium and hydroxyapatite.
- cemented fixation on microblasted surface.
- outer macrostructure: 2 perforated lateral ribs ensuring stability and giving the possibility for a fragments anchoring.
- highly-polished inner surface.
- sterilized by gamma irradiation.

CEMENTLESS RADIAL IMPLANT



Designation

TI/HA RADIAL IMPLANT LS
TI/HA RADIAL IMPLANT LM
TI/HA RADIAL IMPLANT RS
TI/HA RADIAL IMPLANT RM

Ref

MRIHS001
MRIHM001
MRIHS011
MRIHM011

CEMENTED RADIAL IMPLANT



Designation

CEM RADIAL IMPLANT LS
CEM RADIAL IMPLANT RS

Ref

MRICS001
MRICS011

lépine Algérie

9 rue philosophe Tabrizi
Les sources
Bir Mourad Rais - ALGER
lepine-algerie@groupe-lepine.com

lépine Italia

Via Cassanese 100,
20090 Segrate (Milano)
lepine-italia@groupe-lepine.com

lépine Ibérica

J.J. Tadeo Murguía, 3 y 5 Bajos
20304 IRÚN (GUIPÚZCOA)
lepine-iberica@groupe-lepine.com

lépine Maroc

79 avenue IBN SINA
10080 RABAT - AGDAL
lepine-maroc@groupe-lepine.com

175 rue Jacquard - CS 50307
69727 Genay Cedex - FRANCE
Tél. +33 (0)4 72 33 02 95
Fax +33 (0)4 72 35 96 50

www.groupe-lepine.com



groupe lépine

Depuis 1714