

— PROSTHESIS
SCAPHO-TRAPEZIO-TRAPEZOID

INCA™



DOCUMENTATION



groupe lépine

Depuis 1714

INCA™ STT PROTHESIS

- cobalt–chromium–molybdenum alloy.
- cementless fixation: titanium + HA coating.
- anatomical bearing surface with mirror-polished finish.
- keel ensuring primary stability.
- sterilized by Gamma radiation.



STT INCA™

Designation

STT INCA 14/10
STT INCA 16/11

Ref.

MIHC014
MIHC016

INSTRUMENTATION



Designation

STT INCA TEMPLATE


Ref.

MIAGT001



STT INCA IMPACTOR

MIAID001



STT INCA RASP

MIARQ001

lépine Algérie

9 rue philosophe Tabrizi
Les sources
Bir Mourad Rais - ALGER
lepine-algerie@groupe-lepine.com

lépine Italia

Via Cassanese 100,
20090 Segrate (Milano)
lepine-italia@groupe-lepine.com

lépine Ibérica

J.J. Tadeo Murguía, 3 y 5 Bajos
20304 IRÚN (GUIPÚZCOA)
lepine-iberica@groupe-lepine.com

lépine Maroc

79 avenue IBN SINA
10080 RABAT - AGDAL
lepine-maroc@groupe-lepine.com

175 rue Jacquard - CS 50307
69727 Genay Cedex - FRANCE
Tél. +33 (0)4 72 33 02 95
Fax +33 (0)4 72 35 96 50

www.groupe-lepine.com



groupe lépine

Depuis 1714