





BRM[®]
Extremities

Silktoe[®]


Anatomical design
5 right and 5 left sizes
2 pressfit grommets per implant



Right



Left



Hallux rigidus or limitus
Osteoarthritis
Rheumatoid arthritis
Destroyed articular surfaces
Ankylosed joint
Limited range of joint motion



IMPLANT

- **Spacer**
Silicone elastomer
- **Grommets**
Titanium alloy
- **Sterilization**
Ethylene Oxide (EO)

MAIN CONTRAINDICATIONS:

- Active sepsis
- Neurological and/or circulatory and/or skin and/or bone irreparable deficiency and/or non-functioning musculotendinous system
- Psychologically unsuitable patient
- Child patient

CODE

DESCRIPTION

BRDNEWP2020L	Silktoe Spacer - Size 20 Left
BRDNEWP3030L	Silktoe Spacer - Size 30 Left
BRDNEWP3535L	Silktoe Spacer - Size 35 Left
BRDNEWP4040L	Silktoe Spacer - Size 40 Left
BRDNEWP5050L	Silktoe Spacer - Size 50 Left
BRDNEWP2020R	Silktoe Spacer - Size 20 Right
BRDNEWP3030R	Silktoe Spacer - Size 30 Right
BRDNEWP3535R	Silktoe Spacer - Size 35 Right
BRDNEWP4040R	Silktoe Spacer - Size 40 Right
BRDNEWP5050R	Silktoe Spacer - Size 50 Right

INSTRUMENTS

CODE

DESCRIPTION

SET.SILKTOE	Silktoe Instruments Set (full box)
BRDNEWPTSBOX	Silktoe Sterilization Box (no instruments)
BRDNPTHANDLE	Multifunctional Silktoe Handle
BRDNPTHANDSU	Silktoe Handle Sliding Unit
BRDNPTCG2030	20-30 Silktoe Cutting Guide
BRDNPTCG3535	35 Silktoe Cutting Guide
BRDNPTCG4050	40-50 Silktoe Cutting Guide
BRDNPTCGMASK	Silktoe Cutting Guide Mask
BRDNPTR2030P	20-30 Silktoe Metatarsal Rasp
BRDNPTR3535P	35 Silktoe Metatarsal Rasp
BRDNPTR4050P	30-40 Silktoe Metatarsal Rasp
BRDNPTR2030D	20-30 Silktoe Phalangeal Rasp
BRDNPTR3535D	35 Silktoe Phalangeal Rasp
BRDNPTR4050D	40-50 Silktoe Phalangeal Rasp
BRDNPTRADAPT	Silktoe Rasps Adapter
BRDNPTTEST20L	20 Left Silktoe Trial Sizer (Green)
BRDNPTTEST30L	30 Left Silktoe Trial Sizer (Orange)
BRDNPTTEST35L	35 Left Silktoe Trial Sizer (Yellow)
BRDNPTTEST40L	40 Left Silktoe Trial Sizer (Blue)
BRDNPTTEST50L	50 Left Silktoe Trial Sizer (Red)
BRDNPTTEST20R	20 Right Silktoe Trial Sizer (Green)
BRDNPTTEST30R	30 Right Silktoe Trial Sizer (Orange)
BRDNPTTEST35R	35 Right Silktoe Trial Sizer (Yellow)
BRDNPTTEST40R	40 Right Silktoe Trial Sizer (Blue)
BRDNPTTEST50R	50 Right Silktoe Trial Sizer (Red)
BRDNPTI2030P	20-30 Silktoe Metatarsal Grommets Impactor
BRDNPTI3535P	35 Silktoe Metatarsal Grommets Impactor
BRDNPTI4050P	30-40 Silktoe Metatarsal Grommets Impactor
BRDNPTI2030D	20-30 Silktoe Phalangeal Grommets Impactor
BRDNPTI3535D	35 Silktoe Phalangeal Grommets Impactor
BRDNPTI4050D	40-50 Silktoe Phalangeal Grommets Impactor

